KULLANIM KILAVUZU

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T-ENDO MIS NICKEL TITANIUM ROOT CANAL FILE

1. PRODUCT INFORMATION

Only for use in professional dentistry. A dental instrument used for creating an entrance to root canal and expanding and cleaning the space at part of root canal from apical to coronal.

2. PACKING

- T-endo MIS files are offered in sterile packages in blisters of 4.
- T-endo MIS produced with 7 different widths as orificer (25/09), glider (15/04), shaper (20/04), finisher (25/04), (30/04), (35/04) and (40/04)
- Product shelf life is 3 years.

3. INDICATIONS

T-endo MIS files are used for expanding and shaping and mechanical cleaning of root canals during root canal treatment

4. CONTRAINDICATIONS

T-endo MIS files should not be used in root canals with sudden apical curvature that may cause a high risk of fracture.

5. WARNINGS

- T-endo MIS files are sterile. Do not use products with damaged package or opened because their sterilization is lost
- T-endo MIS files are produced only for dentists, suitable for professional use.
- T-endo MIS files cut clockwise, for this reason T-endo MIS files must be used with torque controlled endodontic engines that rotates clockwise move.
- During root canal shaping procedure, grooves on T-endo MIS files must be constantly cleaned
- Do not use files again if you see any deformation on any file.
- T-endo MIS files is suitable to use up to 3 times, be aware to not use it more.
- Ideal speed so that T-endo MIS files can make appropriate root canal expansion is between 250 -300 rpm.
- Ideal torque so that T-endo MIS files can make appropriate root canal expansion is between 2,5-3,5 Ncm.
- T-endo MIS files should be used with light apical pressure; extreme pressure should not be applied towards apical.
- For your own safety, use protective equipment such as gloves, mask and goggles during use of T-endo MIS files.
- The product contains Nickel, this product shouldn't be used for people shows allergy to this material.

6. INSTRUCTIONS FOR APPLICATION

- T-endo MIS "orificer file" Should be used to widen the canal entrance.
- T-endo MIS "glider file" should be used to create a glide path.
- T-endo MIS "shaper ve finisher files" should be used to finish the root canal shaping.
- T-endo MIS files should be used with 2 mm pecking motion, downward and upward
- T-endo MIS files should be made to cut by squeezing into root canals passively
- In case of any deformation during use of T-endo MIS files, they should not continue to be used for root canal shaping procedure, they should be replaced immediately.
- During root canal shaping procedure, when T-endo MIS files are removed from root canal, root canals should be irrigated abundantly.

7. STORAGE CONDITIONS


Files should be stored between 10 and 20 centigrade degrees.

8. STERILIZATION AND DISINFECTION

Files are in sterile blisters. If no deformation is seen after use, files may be used by sterilizing at 134 centigrade degrees for 4 minutes at steam sterilization. Use and sterilization of deformed files with diminished cutting efficiency or used more than three times is not recommended.

All instruments have to be cleaned, disinfected and sterilized prior to each use; this applies to instruments intended to be used again and delivered in sterile condition or opened sterile package. In debt cleaning and disinfection is mandatory for efficient sterilization. In order to clean impurities manually, use a soft, clean brush or soft, clean cloth that would be used only for this purpose; under no circumstances use metal brush or wire wool. Make sure that there is no visible impureness or residue and repeat preliminary cleaning process if needed.

9. SYMBOLS

	Manufacturer
	Date of Production
	Lot Number
	Reference Number
	Consult instructions for use
	Medical device
	Caution
	Use by date
	CE Certificate Notified Body Number
	Keep away from sunlight
	Keep dry
	The product can be sterilized by steam autoclave at 134°C.
	Number of products in the box
	Asorti Blister Package
	Sterilized using irradiation
	Made of Nickel Titanium
	For dental professional use only.
	Do not use if package is damaged and consult instructions for use