PARION FLOW

KOMPOZİT LED IŞIK İLE POLİMERİZE OLAN RADYOOPAK RESTORATİF DİŞ DOLGUSU

1. ÜRÜN BİLGİSİ

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COMPOSITE POLYMERIZED RADIOPAQUE RESTORATIVE DENTAL FILLING WITH LED LIGHT

1. PRODUCT INFORMATION

T-Resto Parion Flow Composite is a self-adhesive, light-curing, resin-based composite that is designed for dental restorative for direct placement. T-Resto Parion Flow Composite eliminates the additional steps of etching/coating/bonding unless otherwise required to bond a resin composition to dentin and enamel. T-RESTO Flow Composite offers advanced bonding effect, high mechanical effect and physical properties similar to other conventional flowable composites.

2. COMPONENTS

T-Resto Parion Flow Composite includes methacrylate resins (Bis-GMA, UDMA, TEGDMA etc.), inorganic filling materials, photoinitiator and inhibitor.

3. PACKAGING

2 pieces of 2 g syringes which are filled with flowable composite and 20 pieces application tips for each syringes are placed in the pouch. Can be preferred in 9 different colors as A1 – A2 – A3 – A3.5 – A4 – B1 – B2 – C1 – C2.

4. INDICATIONS

T-Resto Parion Flow Composite is used for small class I and Class I and Class II base/primer restorations. Additional indications include pitting and cleft tissue adhesives, repair of enamel defects, porcelain restoration, prevention of undercuts, minor occlusal reinforcements in non-stressed areas, and incisors abrasion. It is used in direct restoration of I-II-III-IV-V class cavities, restoration of broken teeth in cavities on the root surface, restoration of diastama and veneer closure cases.

5. CONTRAINDICATIONS

It should not be used in patients with direct pulp capping and sensitivity to any of the product ingredients. In chemically cured composites.

6. WARNINGS AND PRECAUTIONS

- Uncured methacrylate resin can cause contact dermatitis and damage to the pulp. Avoid contact with skin, eyes and soft tissue. In case of contact, rinse thoroughly with water.
- If swallowed, drink plenty of water, do not induce vomiting. Consult a doctor immediately.
- Use personal protective equipment such as disposable masks, gloves and goggles when working with this product.
- Do not use this product for any purpose other than the treatments outlined in the indications section in this user manual.
- This product is intended use for dental purposes only.
- To avoid cross-contamination between patients, the application tips are disposable for each patient.
- Do not use if the syringes are damaged.
- Retract the plunger after application to prevent unwanted flow of paste from the tip.
- It is recommended to use a reliable, high quality, visible light source to cure this product.
- Protect from visible light before use.

7. STORAGE CONDITIONS AND SHELF LIFE

- Syringes should be stored between 15 and 25 degrees Celsius. Should be stored in a dry, cool and dark place. Keep away from direct sunlight.
- The expiry date is indicated on all packaging. Do not use after the expiration date.
- Shelf life is 3 years.

8. APPLICATION

- Determine the desired color in daylight with the VITA scale.
- Teeth must be clean so that color selection is accurate.

- Prepare the cavity for composite restorations.
- Isolate the tooth from the adjacent teeth. Use transparent thin strips.
- The cavity must be clean and dry. Dry the dentin surfaces to strengthen bonding.
- Apply acid and bonding according to the rules.
- Open the cover of the composite east.
- Apply the composite filler not thicker than 2 mm to the dent.
- Keep the beam device as close to the cavity as possible while curing.
- Take care to ensure that the airway is free of oil and other contaminants.
- After the application is finished, close the cap of the tube to make sure that it does not get air.
- Closing the cap of the product after each use does not change the current expiry date. The product can be used during the stated shelf life.
- For repeated use, make sure that the cap is closed correctly after use. After use, the product should be stored in a cool, dry and dark place.

Curing Time;

- 10 seconds for light colors (A1 – A2 etc.)
- 20 seconds for dark colors (A3.5 – A4 etc.)

9. DISPOSED OF

- The disposed products should be disposed of in accordance with hospital/clinical environment specific instructions.

10. SYMBOLS

	Date of manufacture
	Manufacturer
	Consult instructions for use or consult electronic instructions for use.
	Batch number
	Temperature limit
	Caution
	CE Declaration with Notified Body Number
	Keep away from sunlight.
	Keep dry.
	MR Safety
	Medical Device
	Professional for use, only.
	Do not use if package is damaged ad consult instructions for use.
	Product reference code